# Viewing a Device/Devicepool

# Navigation

The Device/Devicepool View page can be accessed via the Devices page or via clicking on a devicepool on the Dashboard page.

## Overview

On this page, device sensor information can be viewed and interacted with.

= Mytridgeonline				ieres - 🔿
0	Techtoard > Denice Paol 3			
🕤 demolijskes so za	Search	AI -		с
E Destinant	064006043004404			
🌲 Alema 🔍 🗸	204800817394134 2 Minuter Age + MK2-681122-0CHONEM			
🗈 Reports 🔍 🗸		1		
🙁 Oneniev 👘 🗸	23.0°C Fridge Temp	1.0°C Room Temp	23.0°C Probe 2	23.0°C Probe 3
🛃 Devices	Normal Min 18 Max 25	Warning Min 18 Max 25	Normal Min 10 Max 25	Normal Min 18 Max 25
🔚 Deviceposis 🗸 🗸	23.0°C Probe 4	23.0°C Probe 5	23 0°C Probe 6	23.0°C Prote 7
🚊 there 🔍 🗸	Normal Vin 15 Max 25	Normal Min 16 Mar 23	Normal Min 18 Max 28	Normal Vie 18 Mar 25
26 Account	- Continue			Training (
	23.0°C Probe 8	80.0% Room Rh	Closed	15.0V Battery
	Normal Min 18 Max 25	Warning Min 42 Max 73	Normal	Warning Mill 2.7 Max 5
	23.0°C Pt100 1	Fault Pt100.2	23.0°C P1100 3	23.0°C Pt100.4
	Normal Win 18 Mix 25	Foult Min 16 Mar 25	Normal Mn 18 Max 28	Normal Vie 18 Mar 25
	00000 01			
	3090P			
	Warning On see Max 1990			
	·			
E+ Logost				

#### **Filtering Sensors**

Sensors can be filtered by their status by using the drop down menu. The search function also allows further filtering.

= Myfridgeonine				<u>tern</u> - Ø
0	n Deshboard > Desice Peol 3			
🕒 demogrikes en an	Search	AI		C
E Deshboard	(	A8		
🌲 Alens 🔍 🗸	264806817394134 3 Minutei Agi + ME2011223-0034004M	Online		
🗈 Reports 🔍 🗸		Offine		
🔿 Ommine 👘 🗸	23.0°C Fridge Temp	1.0 Unset	23.0°C Probe 2	23.0°C Probe 3
📑 Devices	Normal Min 18 Max 25	War-	Normal Min 18 Max 25	Normal Min 18 Max 25
🔚 Deviceposis 🗸 🗸	DD DBD Embed	00.0%0 Proba 5	00.0%0 Broke 6	00.0%0 Pmba7
🛓 there 👘 🗸	23.0 °C	23.0 0	23.0 0	23.0 0
to Account	Normal Min 18 Max 25	Normal Normal Nor 18 Mar 25	Normal Min 18 Max 28	Normal We 18 Mar 25
	23.0°C Probe 8	80.0% Room Rh	Closed Door	15.0V Battery
	Normal Min 18 Max 25	Warning Min 40 Mau 73	Normal	Warning Mill 27 Mail 5
	22 010 P1100 1	Equila Pil012	22 010 Pt100 8	22.0*0 Pt1014
	23.0 0		23.0 0	23.0 0
	Normal win to wat by	Fault 11 Mar 13	Normal 447 10 100 10	Normal on it out of
	3090P Co2			
	Warning Min Son Max 1500			
E+ topost				

### Last Synced Time

The device's last synced time can be viewed here.

= Myfridgeonline				50000 - 🔿
0	Teshboard > Denice Paol 3			
🕤 demolitikes no. n	Search	AI -		c
E Bestboard	264806817204124			
🌲 Alets 🔍 🐱	204800817394134 2 9 Min/H Aya + MK2-681322-00-K0KM			
🗈 Reports 🔍 🐱				
😔 Oreniev 🤍 🗸	23.0°C Fridge Temp	1.0°C Room Temp	23.0°C Probe 2	23.0°C Prote 3
and bevices	Normal Min 18 Max 25	Warning Min 18 Max 25	Normal Min 10 Max 25	Normal Vin 18 Max 25
🔚 Deviceposis 🗸 🗸	23.0°C Probe 4	23 0°C Probe 5	23 0°C Probe 6	23.0°C Prote 7
🚊 Uara 🔍 🗸	Normal Min 18 Mex 25	Normal He 16 Mar 25	Normal Mn 18 Mar 28	Normal Vie 18 Mar 25
Jo Account	realities	Hama	Normal	reating
	23.0°C Probe 8	80.0% Room Rh	Closed	15.0V Battery
	Normal Min 18 Max 25	Warning Min 42 Max 72	Normal	Warning Will 2.7 Max 5
	23.0°C Pt100 1	Fault Pt100.2	23.0°C P1100 3	23.0°C Pt100.4
	Normal Min 18 Mex 25	Foult Min 18 Mar 25	Normal Mn 18 Max 28	Normal Vie 18 Mar 25
	00000 01			
	3090P			
	Warning Un too Max 1300			
Ex Local				

#### **Sensor Information**

Sensors are displayed either on chips or cards. Sensor information can be viewed here such as the sensor name, current value of the sensor, the sensor status and alert configuration. Displayed values are dependent on the template selected. The color of the sensor and the message displayed indicates the sensor's status. This can be one of the following statuses:

• Unconfigured - Grey

- Normal Green
- Min Blue
- LowLow Dark Blue
- Max Red
- HighHigh- Dark Red
- Warning Orange
- Alert Dark Red
- Offline Black
- Fault Dark Green
- Disabled Purple

= Mytridgeonline				tera - Ø.
	Pashboard > Denice Paol 3			
🕤 demo@ikes.co.co	Search	A8		с
E Deshboard	264906817304134			
🌲 Alets 🗠	204600617394134 2 Minute Age + ME2011220CHOHM			
📑 Reports 🗠		9		
🗢 Decriter 🔹	23.0°C Fridge Temp □ Normal Ma 18 Max 25	1.0°C Room Temp Warning Mile 18 Mai 25	23.0°C Probe 2 Normal Mon 10 Max 25	23.0°C Probe 3 Normal Mie 18 Mie 25
E Deviceposis .	23.0°C Probe 4 Normal Min 18 Mar 25	23.0°C Probe 5 Normal He 18 Her 23	23.0°C Probe 6 Normal Mn 18 Max 28	23.0°C Probe 7 Dormal Vie 18 Vie 25
	23.0°C Probe 8 	80.0% Room Rh Warning Min 40 May 70	Closed Door	15.0V Battery Warning Mis 2.3 Max 5
	23.0°C Pt100 1  Normal Min 19 Max 25	Fault Pt100.2 Pault We 18 Mar 28	23.0°C Pt100 3 Normal Mn 18 Max 28	23.0*C Pt100.4 O Normal Win 18 Mar 25
	3090P Co2 Warning Wo co Mex 1500			
E+ Logost				

Alert configurations can be viewed here. The bell icon indications the number of alert groups set, the monitor icon indicates the number of server alerts that are selected and a snooze icon indicates the number of rules that are currently snoozed (this is based on the alert schedule). These are all dependent on the number of rules set on the sensor.

3/8

Last update: 2025/05/01 12:14 user\_manual:dashboard:view\_device https://docs.myfridgeonline.com/user\_manual/dashboard/view\_device

= Mytridgeonine				Serra - 🔿
0	n Deshboard > Denice Peol 3			
🔿 dermo@ikes.co.cs	Search	AI -		С
E Deshboard				
🌲 Alem 👘 🐱	264806817394134 2 Milliona Age + MK2-001322-003400HM			
🗈 Reports 👘 🐱		1		
😔 Overview 🧠 🗸	23.0°C Fridge Temp	1.0°C Room Temp	23.0°C Probe 2	23.0°C Probe 3
🛃 beine	Normal Min 18 Max 25	Warning Min 18 Max 25	Normal Mn 18 Max 25	Normal Min 18 Mai 25
🗄 Deviceposis 🗠 🗸	22 0°C Prote 4	22 0°C Probe 5	22 0°C Probe 6	22.0°C Prote 7
🛓 Uwa 🔍 👒	23.0 C	23.0 C	23.0 C	United 11 No. 25
🛵 Account	Notifial	Norma	NOIMAI	Norma
	23.0°C Probe 8	80.0% Room Rh	Closed	15.0V Battery
	Normal Min 18 Max 25	Warning Min 48 May 79	Normal	Warning Mill 27 Mill 5
	23.0°C Pt100 1	Fault Pt100.2	23.0°C P1100 3	23.0°C Pt100.4
	Normal Win 18 Mex 25	Fault Min 18 Man 15	Normal Min 18 Max 25	Normal Min 18 Mar 25
	3090P			
	Warning Min 500 Max 1500			
E+ topost				

#### **Editing Sensors**

Sensors can also be selected by clicking on them, the sensor color will change to pink to indicate that it has been selected. Upon selection, a side menu will appear with a number of options (depending on the user role).

= ay noted inc	
O Parito and 3	
🔿 demolphes to be	C
St Instituted	
Alets 204800681/394134 Thimman Ap. + MRO 01223 0046948M	
	5
23.0°C Fridge Temp 1.0°C Room Temp 23.0°C Probe 2 23.0°C Probe 23.0°C Probe	#3 ×
Libetes Normal Win 18 Max 25 Warring Win 18 Waa 25 Normal Min 18 Max 25 Normal Win	18 Max 16
Cencepter	x 7
	18 May 25
Normal Norma	
23.0°C Probe 8 80.0% Room Rh Closed Door 15.0V Bern	ery
Normal Min 18 Max 25 Warring Min 46 Min 79 Normal Warring Min	87 Mai 5
23.0°C P1001 Fault P100.2 23.0°C P100.8 23.0°C P10	04
Normal Vin 19 Max 25 Fault Vin 18 Vin 25 Normal Vin 19 Max 25 Normal Vin	18 Max 25
3090P	
Warning We see Mux 1980	
F# Logost	

Selecting Cancel will end the selection.



Selecting Edit will open the Sensor Edit dialog.



Selecting Graph will open the Sensor Graph dialog using the current selected sensor information.



Selecting Delete will delete the selected sensors.



From: https://docs.myfridgeonline.com/ - **MyFridgeOnline** 

Permanent link: https://docs.myfridgeonline.com/user\_manual/dashboard/view\_device

Last update: 2025/05/01 12:14

