


Allocating devices to device pools

Overview

The **Device Pool Management Page** allows users to **view, search, filter, select, and manage devices** within the device pool. From this page, users can assign actions to multiple devices, manage device pools, and easily copy device serial numbers.



Page Overview

When you open the page, you will see two main sections:

- * **Header & Actions Area** – Search, filters, and action buttons for managing devices
- * **Devices Table** – A paginated, sortable table displaying devices in the pool

Searching & Filtering Devices

Basic Search

1. Use the **search bar** at the top of the page.
2. Enter a device-related value (for example, serial number or name).
3. The table refreshes automatically to show matching devices.

Advanced Filters

1. Open the **filter panel** in the header area.
2. Select one or more filters (such as device type or status).
3. Active filters are displayed below the header.
4. The device list updates to reflect your selections.

Selecting Devices

Single & Multiple Selection

* Click the **checkbox** next to a device to select it. * Select multiple devices to perform bulk actions.

Selection Counter

The header displays the **number of selected devices**, helping you track bulk selections.

Clearing Selection

Use the **clear selection** option in the header to deselect all selected devices.

Managing Device Pools

Viewing Device Pools

Each device row displays its associated **device pools** as chips.

* Up to three pools are shown by default. * Additional pools may be hidden if more exist.

Removing a Device from a Pool

1. Locate the device in the table. 2. In the **Device Pools** column, find the pool chip you want to remove. 3. Click the **remove (x)** icon on the chip. 4. Confirm the action if prompted.

The device will be removed from the selected pool.

Copying Serial Numbers

To copy a device serial number:

1. Locate the **Serial Number** column. 2. Click the serial number chip. 3. The serial number is copied to your clipboard. 4. A tooltip or notification confirms the copy action.

Sorting Devices

To sort the device list:

- * Click on a **column header** to sort by that column. * Click again to reverse the sort order.

The table refreshes automatically using the selected sort order.

Pagination

If many devices are available, they are split across multiple pages.

To navigate:

- * Use the **pagination controls** at the bottom of the table. * Select how many rows to display per page (25, 50, or 100).

Loading & Feedback

- * A loading indicator is displayed while device data is being retrieved. * Bulk actions and updates may briefly refresh the table.

Best Practices

- * Use filters to quickly narrow down large device lists. * Double-check selected devices before performing bulk actions. * Use serial number copy functionality to avoid manual typing errors.

Troubleshooting

Issue	Suggested Action
No devices displayed	Clear filters and try searching again.
Unable to remove device pool	Verify you have the required permissions.
Selection not updating	Clear selection and reselect devices.

Related Pages

* Device Details * Device Pool Administration * User Permissions

Last
update:
2026/01/19 07:59 user_manual:devicepools:allocate_devices https://docs.myfridgeonline.com/user_manual/devicepools/allocate_devices?rev=1768809554

From:
<https://docs.myfridgeonline.com/> - **MyFridgeOnline**

Permanent link:
https://docs.myfridgeonline.com/user_manual/devicepools/allocate_devices?rev=1768809554

Last update: **2026/01/19 07:59**

